










	SUPREM LIGHT SHAPED PADS FOR LIGHT TO MODERATE INCONTINENCE					
	Light Mini	Light Normal	Light Extra	Light Extra +	Light Super	
Faecal Containment 870ml to 1000ml						
Very Severe 2800ml plus						
Severe 2150ml to 2800ml						
Heavy 1700ml to 2150ml						
Moderate 1000ml to 1700ml					✓	
Light 340ml to 1000ml		✓	✓			
Stress 200ml to 340ml	✓	✓				
All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes.						
International Standard Method ISO 11948-1						
Code	04020	04207	04208	04210	04034	
CLASSIC FORM SHAPED PADS FOR MODERATE TO HEAVY INCONTINENCE						
	Regular	Extra	Super	Extra +		
Faecal Containment 870ml to 1000ml						
Very Severe 2800ml plus						
Severe 2150ml to 2800ml						
Heavy 1700ml to 2150ml		✓	✓	✓		
Moderate 1000ml to 1700ml	✓					
Light 340ml to 1000ml					✓	
Stress 200ml to 340ml						
All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes.						
International Standard Method ISO 11948-1						
Code	04200	04201	04202	04213		
CLASSIC PAD INSERTS AND PE BACKED PADS FOR REINFORCED ABSORPTION						
	Insert Mini	PE Backed Mini	PE Backed Extra	PE Backed Maxi		
Faecal Containment 870ml to 1000ml						
Very Severe 2800ml plus						
Severe 2150ml to 2800ml						
Heavy 1700ml to 2150ml						
Moderate 1000ml to 1700ml					✓	
Light 340ml to 1000ml	✓	✓	✓			
Stress 200ml to 340ml						
All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes.						
International Standard Method ISO 11948-1						
Code	04024	03991	03993	03995		
SUPREM FOR MEN FOR LIGHT INCONTINENCE						
	Super					
Faecal Containment 870ml to 1000ml						
Very Severe 2800ml plus						
Severe 2150ml to 2800ml						
Heavy 1700ml to 2150ml						
Moderate 1000ml to 1700ml						
Light 340ml to 1000ml					✓	
Stress 200ml to 340ml						
All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes.						
International Standard Method ISO 11948-1						
Code	04034					

	SUPREM PANTS				
	Maxi / S	Extra / M / L / XL			
Faecal Containment 870ml to 1000ml					
Very Severe 2800ml plus					
Severe 2300ml to 2800ml					
Heavy 1700ml to 2300ml	✓				
Moderate 1000ml to 1700ml		✓			
Light 340ml to 1000ml					
Stress 200ml to 340ml					
All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes.					
International Standard Method ISO 11948-1					
Code	04090	04230			
ALL-IN-ONES FOR MODERATE TO SEVERE INCONTINENCE					
	Extra + / X Sm	Extra + / Sm		Regular + / M	Regular + / L
Faecal Containment 870ml to 1000ml					
Very Severe 2800ml plus					
Severe 2300ml to 2800ml				✓	✓
Heavy 1700ml to 2300ml	✓	✓			
Moderate 1000ml to 1700ml					
Light 340ml to 1000ml					
Stress 200ml to 340ml					
All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes.					
International Standard Method ISO 11948-1					
Code	415410SU	04216	04218	04222	
ALL-IN-ONES FOR HEAVY TO VERY SEVERE INCONTINENCE					
	Maxi / Sm	Super + / M		Super + / L	Extra + / L
Faecal Containment 870ml to 1000ml					
Very Severe 2800ml plus		✓		✓	✓
Severe 2300ml to 2800ml					
Heavy 1700ml to 2300ml	✓				
Moderate 1000ml to 1700ml					
Light 340ml to 1000ml					
Stress 200ml to 340ml					
All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes.					
International Standard Method ISO 11948-1					
Code	04217	04220	04224	04226	
ALL-IN-ONES FOR VERY SEVERE INCONTINENCE					
	Max / M / L / XL				
Faecal Containment 870ml to 1000ml		✓			
Very Severe 2800ml plus					
Severe 2300ml to 2800ml					
Heavy 1700ml to 2300ml					
Moderate 1000ml to 1700ml					
Light 340ml to 1000ml					
Stress 200ml to 340ml					
All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes.					
International Standard Method ISO 11948-1					
Code	04221, 04225, 04227				
CLASSIC BED PADS					
	40 x 60cm	60 x 60cm		60 x 75cm	60 x 90cm
Faecal Containment 870ml to 1000ml					
Very Severe 2800ml plus					
Severe 2300ml to 2800ml					
Heavy 1700ml to 2300ml				✓	
Moderate 1000ml to 1700ml			✓		
Light 340ml to 1000ml	✓	✓			
Stress 200ml to 340ml					
All absorbencies quoted are Total Absorbency as determined by a Total Saturation Test over 30 minutes.					
International Standard Method ISO 11948-1					
Code	LFBD8111	04203	04204	04206	